

SUBLIO

HYPERIONIC WATER LIFE

leading
e-biotechnology



6

Whatever has already been invented belongs to the past, only innovations are worthy of interest.

Ettore BUGATTI



This is a revolution.



This is Sublio Hyperionic Water Life.

A unique holistic cosmetic-therapeutic solution to naturally and highly effectively address problems linked to immunosenescence, skin ageing and the majority of skin diseases.



€4 MILLION

4 MILLION EUROS INVESTED



15 YEARS


15 YEARS OF R&D



Hyperionic Water Life technology is a **disruptive** e-biotechnology process based on water hyperionisation.



a world first
revolutionary system



Technology designed
FOR YOUR HEALTH AND WELL-BEING



SUBLIO technology

- does not add anything to the water
- nor does it remove anything from the water




Sublio Hyperionic Water Life is installed in the machine rooms of treatment centres:

- Thermal Baths
- Thalassotherapy Centres
- Balneotherapy Centres
- Spas or Swimming Pools



sublio's actions

Sublio's Hyperionic water



AQUATIC THERAPY PROFESSIONALS

Thermalism - Balneotherapy- Thalassotherapy

Other than the active and visible repair of the skin, Sublio also kicks off an extraordinary action to increase the intracutaneous osmosis of trace elements.

This exceptional process boosts the benefits of aquatic treatments, notably allowing mature skin - often denied the benefits of ionic water due to “short” dermal vein circuits - to be able to capture the benefits of trace elements - whose highly beneficial effects on the immune system are well-known - and to circulate them throughout the body.



scientifically
proven
effectiveness

In order to gauge the efficacy of treatments offered by Sublio hyperionic water, scientific studies have been conducted by two leading laboratories.

The results of ex vivo objectification studies combined with genome studies through modulated messenger RNA transcriptome analysis demonstrate Sublio water's **unrivalled efficacy** in the anti-ageing process by actively repairing the skin barrier, as well as maximising penetration by trace elements and cosmetic-therapeutic active ingredients.

For its biotechnology certifications, Sublio France has chosen two leading laboratories : BIO-EC & GENEX.

LABORATOIRE
BIO-EC
Centre de recherches biologiques
& d'expérimentations cutanées

GENEX
GENOMIC SOLUTIONS



in detail

deep repair and regeneration of the skin



(...) SUBLIO Hyperionised water [editor's note: although applied in very small quantities (2mg/cm²)] completely stopped skin drying out from the first application.


(...) Sublio Hyperionised water has a clear moisturising effect.

(...) A single application of SUBLIO water offers the same repairing power as a complex cosmetic formulation.


(...) The hyperionic water produced using the SUBLIO HYPERIONIC WATER LIFE system has an excellent repairing effect on the skin.

(excerpts from the study published by Laurent Penno-Mazzarino, Technical Director at the Bio-EC Laboratory)

Schematic representation of the skin's condition




without Sublio water




with Sublio water


In-depth histological analyses, based on changes to marker protein activity in epidermis reconstruction (involucrin, filaggrin, aquaporin 3) were conducted to highlight the hydration and reconstruction of the skin's protective layer (stratum corneum).




SUBLIO's hyperionic water causes the skin's outer layer to become thicker by increasing layering.



After skin delipidation, SUBLIO water prevents a decline in filaggrin levels by doubling its initial rate after 24 hours.



SUBLIO hyperionised water offers clear hydrating activity confirmed by the reduced expression of aquaporin 3.



The clear overexpression of involucrin suggests that SUBLIO water actively plays a role in reinforcing the stratum corneum.




in detail

maximizing penetration by trace elements



Y O U N G . S K I N



Trace element penetration without Sublio


The skin barrier is altered, the strength of the transcutaneous electrical current is low, intradermal penetration of trace elements is moderate and shallow.

Trace element penetration with Sublio

The skin barrier is totally repaired, the strength of the transcutaneous electrical current is powerful, the intradermal penetration of trace elements is significant and deep.



M A T U R E . S K I N



Trace element penetration without Sublio

The skin barrier is altered, the strength of the transcutaneous electrical current is low, the intradermal penetration of trace elements is moderate and shallow without being able to reach the "short" dermal vein circuits.

Trace element penetration with Sublio

The skin barrier is totally repaired, the strength of the transcutaneous electrical current is powerful, the intradermal penetration of trace elements is significant and sufficiently deep to reach the "short" dermal vein circuits.



| the sublio france experts



Frédéric ESNAULT

Chief Executive Officer
& Founder



Abdellatif BAABOUZ

Chief Technical Officer
& Co-Founder



Pascal LEBLOND

Chief Product Officer
& Co-Founder



Laurent MARECHAL

Chief Operating Officer
& Co-Founder

| a message for you



“

Having read this brochure, you can probably appreciate the extensive benefits offered by our breakthrough technology, whose scientifically proven effectiveness establishes a totally new era in the therapeutic performance of hydrotherapy.

Frédéric Esnault

Frédéric ESNAULT
CEO, Sublio France



SUBLIO

HYPERIONIC WATER LIFE PRO

NON-CONTRACTUAL ILLUSTRATION



SUBLIO

HYPERIONIC WATER LIFE

NON-CONTRACTUAL ILLUSTRATION



QUICK AND EASY TO INSTALL

1

Recommendation

Depending on your hyperionic needs, we determine the positions and assess the number of machines required to implement your project.

2

Installation

Our fitters transport, position and start up Sublio with the support of your professionals for the required electrical and hydraulic connections.

2 days

Depending on the specific features unique to the installation of each Hyperionic Water Life, the positioning and connections generally require less than two days of work on site.

→ Our expert partners

PureFlow®

PLASTI-PLAN
B.E MOULES ET PIÈCES PLASTIQUES

RAYON
USINAGE
Mécanique de précision

SAT

dv ■ Design
■ Inox
■ Vendée

OUR GUARANTEES

1

The scientific results are 100% certified.

2

The Hyperionic Water Life generator has a **12-year parts and labour warranty** excluding flood, fire or power surge.

3

The use of this **warranty** requires the replacement of the device's two induction plates as part of its **annual maintenance**.

4

The Sublio range is entirely made in France.

5

Sublio customer services are available **24/7**.



SUBLIO

HYPERIONIC WATER LIFE



SUBLIO FRANCE
8 rue René Coty - CS 40 747
85 018 LA-ROCHE-SUR-YON CEDEX
TEL: +33 (0) 9 72 10 15 15

sublio.com
hello@sublio.com